



A Sensitive and Selective GC-MS Method for Analysis of Genotoxic Impurities in Dobutamine Hydrochloride

VIJAYA LAKSHMI.MADDALA*¹, P.C.RAY¹ and KMVN RAO¹

¹Inogent Laboratories Private limited.

*Corresponding author E-mail: Vijaya_chem@yahoo.com

<http://dx.doi.org/10.13005/ojc/320347>

(Received: April 06, 2016; Accepted: May 27, 2016)

ABSTRACT

A sensitive GC-MS method for analysis of residues of Methyl-p-toluene sulfonate [Fig.2], Ethyl p-toluene sulfonate [Fig.3] and Isopropyl p-toluene sulfonate [Fig.4] Genotoxic impurities in Dobutamine hydrochloride [Fig.1] drug substance, has been developed and validated. Under the optimum conditions of the proposed method is specific, precise, linear, accurate and rugged. The response was linear in the concentration range of LOQ to 0.18 mg.g⁻¹ with correlation coefficients of greater than 0.99 for each impurity. Under the optimum conditions used for the sample preparation the recovery of each impurity is in between 78-108%.

Keywords: Dobutamine hydrochloride, GC-MS, Genotoxic impurities, method validation.

INTRODUCTION

Genotoxic impurities have the tendency to cause a potential damage of the genetic material, for example DNA even at very low level of exposure. There are three primary effects that genotoxins can have on organisms by effecting their genetic information. Genotoxins can be carcinogens, or cancer causing agents, mutagens or mutations-causing agents or teratogens, birth defect causing agents. Thus, there is a high demand and need to identify the potential impurities that arise during synthesis, purification and storage keeping the damaging effects of Genotoxic impurities, many

international regulatory agencies has already published the guidelines regarding the limits of Genotoxic impurities. Both EMEA guidelines and a phRMA propose a maximum of daily exposure target of 1.5µg per day genotoxic impurities in pharmaceuticals is recommended Based on the Threshold of toxicological concern (TTC), the concentration limits of the impurity in a drug substances or drug products can be derived based on the maximum daily dose: concentration (ppm) = [1.5µg/ day]/[dose (g/day)]. For a drug dosed at 1g per day, for example, 1.5 ppm would be the limit of that specific Genotoxic impurity.

Dobutamine hydrochloride¹⁻² is a sympathomimetic with direct effects on beta1- adrenergic receptors, giving it a prominent isotropic action on the heart. It also has some alpha – and beta2 – agonist properties. Although it is structurally related to dopamine it has no specific dopaminergic properties; however, like dopamine, the isotropic action of dobutamine on the heart is associated with less cardiac – accelerating effect than that of isoprenaline. Dobutamine is used to increase the contractility of the heart in acute heart failure, as occurs in cardiogenic shock and myocardial infarction. It is also used in septic shock. Other circumstances in which its isotropic activity may be useful are during cardiac surgery and positive end – expiratory pressure ventilation.

During the synthesis these are the three Genotoxic impurities present in the Dobutamine hydrochloride. Based on the literature survey, to the best of our knowledge there is no method to estimate of these Genotoxic impurities. Here in this paper we present a sensitive GC-MS³ method has been developed and validated for the estimation of Genotoxic impurities. The mentioned Genotoxic impurities are prepared and tested in Inogent Laboratories private limited, Nacharam, Hyderabad.






Fig. 1: Dobutamine hydrochloride



Chemical Formula: C₈H₁₀O₃S
Molecular Weight: 186.23

Fig. 2: Methyl p-toluene sulfonate



Chemical Formula: C₉H₁₂O₃S
Molecular Weight: 200.25

Fig. 3: Ethyl p-toluene sulfonate

EXPERIMENTAL

Active pharmaceutical ingredient standards and samples were supplied by Inogent Laboratories private Limited, Hyderabad, India. The HPLC grade Toluene was purchased from Merck, Darmstadt, Germany.


Chromatographic conditions and equipment

Mutagenic impurities⁵⁻⁶ are separated by GCMS on HP-5 ms 30mX0.25mm, 0.25μ column, Samples and standards are prepared in 0.1% Diethyl amine in methanol in the system Agilent 7890B Gas chromatograph equipped with 5977A Mass selective detector. The oven programme consists of Initial the oven is programmed at 60°C for 3 minutes, then increased at the rate of 25°C/min upto 250°C, hold the temperature at 250°C for 10.5 minutes, Injector temperatue is put at 200°C, Auxiliary temperature at 300°C. The carrier gas is Helium(High purity 99.999%), at a constant flow rate of 1.4 mL/min and the split ratio is 1:10. having the injection volume 1.0 μL, the total runtime of 21 minutes.

The Ion source in the mass spectrophotometer is the Electron impact ionization, in which the acquisition mode is SIM(Selective ion monitoring), having the source temperature at 230°C and quadrapole temperature at 150°C, solvent delay time is put for 7 minutes. The SIM parameters are Group id is 1, high resolution and the ions/dwell time are 155/100, 172/100, 186/100, 200/100, 214/100.

Preparation of the stock solutions

An individual stock solution (1000 $\mu\text{g}\cdot\text{mL}^{-1}$) of all impurities was prepared in diluent. Sample



Chemical Formula: $\text{C}_{10}\text{H}_{14}\text{O}_3\text{S}$
Molecular Weight: 214.28

Fig. 4: Isopropyl p-toluene sulfonate



solution was prepared ($10 \text{ mg}\cdot\text{mL}^{-1}$) in diluent. Limit of detection and limit of quantification for the impurities were evaluated by preparing the impurities in the level from $60 \mu\text{g mL}^{-1}$ to $150 \mu\text{g mL}^{-1}$. Linearity of each impurity is determined from LOQ level to $180 \mu\text{g mL}^{-1}$.

RESULTS AND DISCUSSION

Method development

Method development was planned for the estimation of the above mentioned Genotoxic impurities in Dobutamine hydrochloride by using Gas chromatographic technique with mass spectrophotometer as detector, as evaluation limit was as low as $120 \mu\text{g mL}^{-1}$. Different columns with different stationary phases and dimensions were used like DB-624, DB-1, RTX-624 and RTX-1301 for

System suitability chromatogram



the analysis. Finally separation is achieved on HP-5 MS column, mid polar column with 5% phenyl and 95% dimethyl polysiloxane as a stationary phase. During the method development using this column, appropriate separation is achieved between the impurities.

Method validation

Method validation¹¹ activity was conducted with necessary parameters as per the analytical method validation guidelines of International conference on Harmonization.

Blank chromatogram



Sample chromatogram




Table 1: LOD/LOQ

LOD/LOQ	Methyl p- toluene sulfonic acid	Ethyl p- toluene sulfonic acid	Isopropyl p- toluene sulfonic acid
Limit of quantification (ppm)	13	8	12
Limit of detection(ppm)	41	24	36

Table 2: Precision at limit of quantification level

Parameter	Methyl p- toluene sulfonic acid	Ethyl p- toluene sulfonic acid	Isopropyl p- toluene sulfonic acid
%RSD for six preparations	3.3	5.9	5.7

Table 3: Linearity

Parameter	Methyl p- toluene sulfonic acid	Ethyl p- toluene sulfonic acid	Isopropyl p- toluene sulfonic acid
Correlation coefficient	0.997	0.997	0.997
Slope	262.993	612.167	558.855
Y-Intercept	2746.472	-525.528	-6081.424

Table 4: Accuracy

Accuracy Level	%Recovery		
	Methyl p- toluene	Ethyl p- toluene	Isopropyl p- toluene
LOQ Level	96.9	84.7	89.1
50% level(60 ppm)	96.0	94.4	90.7
100% level(120 ppm)	83.9	91.4	89.8
150% level(180 ppm)	87.2	78.3	92.3

Table 5: Method precision and Ruggedness

Parameter	%Recovery		
	Methyl p- toluene	Ethyl p- toluene	Isopropyl p- toluene
Method precision %RSD for six preparations	4.5	2.7	3.2
Ruggedness % RSD for six preparations	5.7	6.9	7.1

Specificity

Specificity is done satisfactorily by checking the all impurities individually and with spiked solution, confirmed that there is no interference of any impurity.

Limit of detection (LOD) and Limit of Quantification (LOQ)

The LOD and LOQ values for each impurity were predicted by using the series of solutions from 60 $\mu\text{g mL}^{-1}$ to 150 $\mu\text{g mL}^{-1}$ by using standard deviation of the response and slope method. The limit of detection and limit of quantification results are as shown below table-1. The method is precise at Limit of quantification level. The respective results are represented in below table-2

Linearity⁹

Linearity of the three impurities were satisfactorily done. A series of solutions were prepared at test concentration levels from around LOQ level to 180 $\mu\text{g mL}^{-1}$ level. The peak area versus concentration data was done by linearity plot slope, intercept and residual sum of squares analysis. The calibration curve given based on response over the concentration range of the impurities. The results are tabulated in table-3

Precision

The precision⁸ of the developed method was checked by preparing the solutions by spiking the impurities at 100% level(120 $\mu\text{g mL}^{-1}$) with the drug substance for six times to show that the system is precise. The %Relative standard deviation of the

areas at each level is less than 8% confirming the good precision. The method is also proved as rugged by analyzing the same with different column and different analyst on different day. The results are represented in table-5

Accuracy

The accuracy of the method was evaluated in sample solutions were prepared in triplicate by spiking all impurities at LOQ level, 60 $\mu\text{g mL}^{-1}$, 120 $\mu\text{g mL}^{-1}$ and 150 $\mu\text{g mL}^{-1}$. The % recoveries of the impurities are in table-4.

ACKNOWLEDGMENTS

The author grateful to Inogen Laboratories private Limited for providing the necessary facilities during the course of investigation.

CONCLUSION

In this paper a sensitive specific, accurate, validated and well defined GC-MS method for the quantification of Genotoxic impurities in Dobutamine hydrochloride at ppm level was described. The described method is highly reliable technique for the quantification of impurities present in the Dobutamine hydrochloride.

REFERENCES

1. Thippani ramesh, Nageswararao pothuraju, Ramisetty Nageswararao, saida shaik, *Journal of pharmaceutical analysis*, **2013**, *3* 434-439.
2. Kondamadgu S, Gandu V. Quantification of Drugs and Pharmaceuticals Using N-Bromosuccinimide and Methyl. *Orient J Chem* **2016**;32(1).
3. Sahil K, Prashant B, Akanksha M, Premjeet S, Devashish R ()GC-MS: Applications. *International Journal Pharma & Biological Archives* **2011** *2*: 1544-1560.
4. Senthamil Selvan P, Veeran Gowda K, Mandal U, SamSolomon WD, Pal TK, *J Chromatogr B*, **2007**, *858*, .,143-150.
5. Kamila MM, Mondal N, Ghosh LK, Gupta BK, *Pharmazie*, **2007**, *62*,486-487.
6. Maurer HH, Tenberken O, Kratzch C, Weber AA, Peters FT *J Chromatogr A*, **2004** , 169-181.
7. Ramakrishna NVS, Vishwottam KN, Koteswara MS, Santosh MM, Varma DP, *J Pharm Biomed Anal* , **2005**, *39*, 1006-1013.
8. <http://www.emea.eu.int/pdfs/human/swp/519902en.pdf> .
9. <http://www.fda.gov/downloads/Drugs/GuidanceComplianceRegulatoryInformation/Guidances/ucm079235.pdf>.
10. M. Sun, L. Bai, and D. Q. Liu, *Journal of Pharmaceutical & Biomedical Analysis*, **2009**.49(2): p.529-533.
11. International Conference on Harmonization of Technical Requirements for Registration of Pharmaceuticals for Human Use (I.C.H.), Q2/R1, Validation of analytical procedures Text and methodology, 1995.